



PRIDE[®]
MOBILITY
LIVE YOUR BEST[®]

Tel: +31 (0)85-0408630 . Fax: +31 (0)85-0430039
pridemobility.eu . info@prideeurope.nl


De Victory Luna scootmobiel met 3 wielen levert krachtige prestaties zowel binnen als buiten.



13" luchtbanden



Eenvoudig te bedienen display



Voormand standaard

Specificaties

Maximale belasting	159 kg
Snelheid	Tot 15 km/u
Maximale actieradius	Tot 41 km
Draaicirkel	1257 mm
Bodemvrijheid	95 mm
Lengte	1397 mm
Breedte	699 mm
Banden	13" luchtbanden, voor en achter
Aandrijving	Achterwielaandrijving, afgedichte transaxle, 24 volt gelijkstroommotor
Gewicht zonder accu's	87,5 kg
Accu's	Type: Twee 12-volt, deep-cycle Capaciteit: 55 Ah of 75 Ah
Meer?	Download het specificaties blad op de website

Kenmerken

- Gemakkelijke gasbediening door schakelbare duw- en trekhendel
- Grote 13" luchtbanden met grip op het terrein voor uitstekende prestaties buitenshuis
- Volledig verlichtingspakket inclusief richtingaanwijzers en in hoek verstelbare LED-kop- en achterlichten
- Voormand standaard
- Helmstokdisplay met slechts één hand te bedienen

1. Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition, and tire condition. This specification can be subject to a variance of (+ or -) 10%.
2. Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or -) 3%.
3. Battery weight may vary based on battery manufacturer.
4. Includes standard seat weight.

NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%. Dimensions listed are for power base only. Overall measurements will vary based on seating and accessory selections.

Theoretical range is calculated under test conditions in accordance with European standards. Actual range in normal use will depend on many factors, including the condition of the vehicle and its batteries, the weight of the driver, correct tyre pressures, ambient temperature, and the gradient and surface of the road or pavement. Dimensions listed are for power base only. Overall measurements will vary based on seating and accessory selections.

Beschikbare kleuren



Rood